



REGENYAL idea

Safe treatment for face re-shaping


REGENYAL IDEA is a re-absorbable bio-revolumetry product (sterile, apyrogenic and physiological gel) used as filling material in corrections of medium kind of skin sagging of the face. The main ingredient is Hyaluronic Acid obtained through bacterial fermentation. **Regenyal Idea** spreads naturally within the areas surrounding the implantation, offering this way a balanced and continuous distribution of the substance.

The Hyaluronic Acid implemented is viscous and elastic to a degree that enables the volume to remain compact without interfering with the ability of the latter to adapt to individual morphology and to maintain natural features. Thanks to its special properties, such as **hydration, diffusivity, physiological correction** and **lifting effect**, Regenyal Idea is indicated for those who wish to:

- RESTORE VOLUMES
- CORRECT THE LOWER ZYGOMATIC AREA
- CORRECT THE NASOGENIAN AREA
- CORRECT THE MOUTH ANGLES
- CORRECT THE CHIN AREA
- REDUCE MEDIUM WRINKLES



REGENYAL
LABORATORIES



BIOREVOLUMETRY

Manufacturing Technology

REGENYAL IDEA is a mono-phase gel obtained through a long-lasting cross-linking process at low temperature (gentle cross-linking). Slow mixing enables the cross-linking agent (BDDE) to spread better and in a more uniform manner.

The amount of BDDE used is 30 per cent lower than in the majority of other products, available on the market, having the same degree of viscosity.

In addition, two molecular weights of Hyaluronic Acid are used, at 1M and at 2M Dalton, which allows for the best calibration of the product's diffusivity properties.

Highly Pure Hyaluronic Acid

The Hyaluronic Acid used for the manufacturing guarantees for a maximum quality standard in terms of raw material purity; in particular, the complete absence of nucleic acids and proteins stands for **SUPERIOR QUALITY** and **ABSOLUTE SAFETY**.

The degree of sophistication obtained in our laboratories enables us to obtain a highly efficient and reliable bio-revolumetry product.

PYRAMID-SHAPED IMPLANTATION - Zygomatic area correction




Ingredients

- Cross-linked hyaluronic acid.....25 mg/g
- Phosphate buffer;
Water for injection solutions.....1g

The innovative manufacturing method is covered by an **INTERNATIONAL PATENT**

REGENYAL IDEA - RECOMMENDED IN CASE OF LOCAL HYPO-VOLUMES (LH)

Mode of use



It is recommended to repeat the treatment after 7/8 months

NASOGENIAN WRINKLES - MOUTH ANGLES

